

# Temp-Secur® Collection and Transport Bottles

These 100 ml Urine Specimen Collection Bottles feature a patented design that helps prevent tampering by deterring the use of substituted or externally heated specimens. Each bottle includes a concealed temperature monitoring thermometer hidden beneath a tear-strip label, allowing collectors to verify specimen temperature only after collection.



FS005

## SPECIFICATIONS

### Construction

- Concealed temperature monitoring thermometer
- Heat-sealed integrity pouch
- External label w/ tear-strip

### Performance

- Transport and collect urine samples
- Prevents subjects from using someone else's urine

### Capacity

- 100 ml

### Quantity

- 100 per case

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI