

SPECIMEN COLLECTION

120 mL Urine Specimen Collection and Transport Jar

Our 120 mL pre-sealed collection jar features a wide 2-inch opening, making it especially convenient for female users. Urine Specimen Collection and Transport Jars are available with the option of an attached temperature strip or without one.



FS006 and FS007

SPECIFICATIONS

Construction

- Wide 2" opening
- Optional temperature strip for quick verification

Performance

- Transport and collect urine samples
- Ideal for female subjects

Capacity

- 120 ml

Quantity

- 150 per case

Variations

Model No.	Description
FS006	Without Strip
FS007	With Strip

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI