

240 mL Urine Specimen Collection Cup

These urine specimen cups feature a leak-resistant snap-lock lid and are individually packaged in heat-sealed integrity pouches. Each 240 ml cup is available with or without an external urine temperature monitoring strip attached.



FS008 and FS009

SPECIFICATIONS

Construction

- Snap and lock lid
- Leak-resistant
- Wide base and mouth design
- Comes in a heat-sealed integrity pouch
- Available with or without a urine temperature monitoring thermometer affixed to the outside of the cup

Capacity

- 240 ml

Dimensions

- Base: 2.5" diameter
- Mouth: 3" diameter

Quantity

- 140 per case

Performance

- Ideal for on-site testing
- Facilitates urine specimen collection from female subjects

Variations

Model No.	Description
FS008	Without Strip
FS009	With Strip

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI