

Liquid Absorbing Sheets

These liquid absorbing sheets feature an immobilizing agent that reacts on contact with blood, urine, or other body fluids. Each sheet can absorb anywhere from a few drops up to 200 ml of liquid, far exceeding the capacity of standard liquid absorbing packets.



FS010

SPECIFICATIONS

Construction

- Liquid absorbing sheet

Dimensions

- 4.25" x 2.25"

Capacity

- 200 ml

Quantity

- 100 per pack

Performance

- Immobilizing agent gels body fluid on contact
- Absorbs up to 200 ml of liquid
- Collects blood and urine

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI