

# Small Ozone Generating - Air Purifiers

These Ozone Generating Air Purifiers create a controlled, safe level of ozone by converting oxygen molecules into ozone. The extra oxygen atom from ozone attaches to odor-causing molecules, eliminating odors and fumes through oxidation, and leaving the air smelling fresh.



FS011 and FS012

## SPECIFICATIONS

### Construction

- Baked gray enamel metal cabinet
- 57 CFM (cubic feet per minute) fan

### Dimension

- 10.5" L x 6.5" W x 8" H

### Options & Accessories

- FS017: Disposable Sporax Antimicrobial Replacement Filter

### Performance

- Adjustable ozone output control provides from 3 to 300 mg/hr
- Used in small to medium size laboratories
- Up to 10,000 cubic feet
- Eliminates various odors

### Note:

- 5-year warranty included

### Variations

| Model No. | Utilities/ Power |
|-----------|------------------|
| FS011     | 120V/60Hz        |
| FS012     | 230V/50Hz        |

**MOPEC**  
GROUP

**Mopec**  
ELEVATING ANTIBIOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI