

# Medium Ozone Generating - Air Purifiers

These Ozone Generating Air Purifiers create a controlled, safe level of ozone by converting oxygen molecules into ozone. The extra oxygen atom from ozone attaches to odor-causing molecules, eliminating odors and fumes through oxidation, and leaving the air smelling fresh.



FS013 and FS014

## SPECIFICATIONS

### Construction

- Baked gray enamel metal cabinet
- 102 CFM (cubic feet per minute) fan
- Multiple option timer

### Dimension

- 11.75" L x 7.5" W x 8.75" H

### Options & Accessories

- FS017: Disposable Sporax Antimicrobial Replacement Filter

### Performance

- Adjustable ozone output control provides from 190 to 625 mg/hr
- Used in medium size laboratories Up to 40,000 cubic feet
- Eliminates various odors

### Note:

- 5-year warranty included

### Variations

Model No.	Voltage
FS013	120V/60Hz
FS014	230V/50Hz

**MOPEC**  
GROUP

**Mopec**  
CLEANING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI