

Medium Ozone Generating - Air Purifiers

These Ozone Generating Air Purifiers create a controlled, safe level of ozone by converting oxygen molecules into ozone. The extra oxygen atom from ozone attaches to odor-causing molecules, eliminating odors and fumes through oxidation, and leaving the air smelling fresh.



FS013 and FS014

SPECIFICATIONS

Construction

- Baked gray enamel metal cabinet
- 102 CFM (cubic feet per minute) fan
- Multiple option timer

Dimension

- 11.75" L x 7.5" W x 8.75" H

Options & Accessories

- Disposable Sporax antimicrobial replacement filter (FS017)

Performance

- Adjustable ozone output control provides from 190 to 625 mg/hr
- Used in medium size laboratories Up to 40,000 cubic feet
- Eliminates various odors

Note:

- 5 year warranty included

Variations

Model No.	Voltage
FS013	120V / 60Hz
FS014	230V / 50Hz

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI