

# Large Ozone Generating - Air Purifiers

These Ozone Generating Air Purifiers create a controlled, safe level of ozone by converting oxygen molecules into ozone. The extra oxygen atom from ozone attaches to odor-causing molecules, eliminating odors and fumes through oxidation, and leaving the air smelling fresh.



FS015 and FS016

## SPECIFICATIONS

### Construction

- Baked gray enamel metal cabinet
- 102 CFM (cubic feet per minute) fan
- Multiple option timer

### Dimension

- 11.75" L x 7.5" W x 8.75" H

### Options & Accessories

- Disposable Sporax antimicrobial replacement filter (FS017)

### Performance

- Adjustable ozone output control provides from 340 to 1250 mg/hr
- Used in large laboratories  
Up to 60,000 cubic feet
- Eliminates various odors

### Note:

- 5 year warranty included

### Variations

| Model No. | Voltage     |
|-----------|-------------|
| FS015     | 120V / 60Hz |
| FS016     | 230V / 50Hz |

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI