

30" Reveal Fuming Chamber

The Reveal Ductless Automatic Cyanoacrylate Fuming Chamber is a professional forensic tool for safe, consistent fingerprint development. Its durable polypropylene design, Printbuster™ Pro LCD controls, and Multiplex™ Filtration System ensure precise fuming while protecting users from harmful fumes, with secure, customizable features for reliable evidence processing.



FS034 and FS037

SPECIFICATIONS

Construction

- Includes a filtration system:
 - (1) FS040 electrostatic pre-filter
 - FS034: (1) FS035 carbon filter
 - FS037: (1) FS038 carbon filter
- 2 Hanging rods
- Compact fluorescent lighting
- Sturdy polypropylene

Materials

- Tempered Glass
- Ebmpapst™ centrifugal fan

Dimensions

- FS034:
 - External dimensions: 30" W x 48" H x 28" D
 - Internal dimensions: 36" H
- FS037:
 - External dimensions: 30" W x 74.88" H x 28" D
 - Internal dimensions: 57.24" H

Utilities / Power

- 120V / 60Hz or 220V / 50Hz

Options & Accessories

- FS034: FS036 Mobile Cart
- FS037: FS039 Mobile Cart

Performance

- Manage full fuming cycle: humidity, temperature, circulation fan, hot plate, and purge
- Customizable sequences for consistent results
- Noise/dBA/1 meter - < 53
- 145 CFM airflow, 80 FPM face momentum

Variations

Model No.	Description
FS034	30" Reveal Benchtop Fuming Chamber
FS037	30" Reveal Freestanding Fuming Chamber

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI