

Carbon Filter for 30" Fuming Chambers

The carbon filter is designed to capture and neutralize airborne chemical fumes generated during cyanoacrylate fuming. It is specially treated and engineered with a high surface area and controlled porosity to absorb a range of aerosolized chemicals, helping prevent vapors from being released back into the lab environment.



For reference only: Reveal Fuming Chamber

SPECIFICATIONS

Construction

- (1) filter provided upon purchase of the reveal fuming chamber
- Carbon filter used inside the reveal fuming chambers

Utilities / Power

- 2 options:
120V/60Hz
220V/50Hz

Performance

- Check every 30 days
- Captures and neutralizes a range of airborne chemical fumes

Variations

Model No.	Description	Type
FS035	Filter for 30" Fuming Chamber	Benchtop
FS038	Filter for 30" Fuming Chamber	Freestanding

MOPEC
GROUP

Mopec
ELEVATING ANTHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI