


Electrostatic Pre-filters for 30" Fuming Chambers

A pre-filter designed for the Reveal Fuming Chamber that captures airborne particles and aerosols before they reach the main filtration system. It helps extend the life of downstream filters while maintaining efficient airflow and safe workstation operation.



FS039

For reference only: Reveal Fuming Chamber

SPECIFICATIONS

Construction

- (1) filter provided upon purchase of the reveal fuming chamber
- Electrostatic pre-filter used inside the reveal fuming chamber

Utilities / Power

- 2 options:
120V/60Hz
220V/50Hz

Options & Accessories

- FS034: 30" Benchtop fuming chamber
- FS037: 30" Freestanding fuming chamber

Quantity

- 6 per package

Performance

- Check every 30 days
- First layer of filters that uses static electricity to trap particles in the air before reaching the main HEPA filter
- Helps the workstation last longer and perform more efficiently

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI