

# 36" Reveal Ductless Fume Hood

The ductless fume hood provides a safe, efficient workspace with excellent visibility and user-friendly controls. Its advanced filtration system and thoughtful design protect operators while allowing flexible, plug-and-use installation.



FS040

## SPECIFICATIONS

### Construction

- (2) 15 Watt bulbs
- Ebmpapst centrifugal fan
- (1) electrostatic filter
- (1) HEPA filter
- (1) carbon filter
- Manual speed controller
- Hour counter
- Electronic blower switch

### Utilities / Power

- 120V / 60Hz or 220V / 50Hz

### Performance

- Noise / BA / 1 meter - < 53
- Airflow CFM - 206; Face Momentum FPM - 100
- Clear back and side panels as well as a double-hinged self-locking clear front sash

### Dimensions

- External Dimensions: 36" W x 27" H x 31" D
- Internal Dimensions: 19" H

### Materials

- Epoxy coated steel or polypropylene

### Options & Accessories

- FS041: carbon filter
- FS042: HEPA filter
- FS043: Electrostatic pre-filters

### Compatibility

- FS044: Mobile cart

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI