

Medium-Duty Body Bags

These medium-duty body bags are made from layered polyethylene with reinforced scrim for added strength and durability, while remaining flexible in extreme temperatures. Their chlorine-free construction ensures reliable containment and eliminates the disposal concerns associated with vinyl materials.



FS070

SPECIFICATIONS

Construction

- Heavy duty
- (2) layers of 3 mm polyethylene film laminated to 1000 denier scrim reinforcement

Materials

- 3 mm chlorine-free polyethylene
 - 1000 denier scrim reinforcement
- Equivalent strength to 12 mm vinyl

Performance

- Does not release chlorine gas allowing for easier disposal compared to vinyl material
- Stays flexible up to - 70° F

Variations

Model No.	Dimensions	Quantity
FS070	22" X 30"	12 per case
FS071	36" X 60"	12 per case
FS072	36" X 96"	12 per case
FS073	48" X 100"	6 per case

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI