

Single-Pouch Specimen Bag

Tamper-evident specimen bags designed to safely contain leaks from urine bottles. Each bag features a secure adhesive seal that cannot be reopened without damage and is available with optional absorbent sheets that gel on contact with liquid, effectively handling anything from small spills to full bottle contents.



SPECIFICATIONS

Construction

- Tamper-evident seal
- Adhesive strip

Dimensions

- 6" L X 6" W

Performance

- Contains liquid in case of urine bottle leak
- Sheets absorb up to an entire contents of specimen bottle

Quantity

- 200 per package

Variations

Model No.	Description
FS081	Without Absorbing Sheet
FS082	With Absorbing Sheet

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI