

Double-Pouch Specimen Bag

Specimen bags designed to securely contain leaks from urine bottles, featuring a tamper-evident adhesive seal that cannot be reopened without damage. With dual compartments to separate specimens from paperwork—and optional absorbent sheets that gel on contact with liquid, handling anything from minor spills to full bottle contents.



SPECIFICATIONS

Construction

- Tamper-evident seal
- Adhesive strip
- (2) Compartments

Dimensions

- 6" L X 10" W

Performance

- Contains liquid in case of urine bottle leak
- Sheets absorb up to an entire contents of specimen bottle
- Separates specimen bottle(s) from documentation

Quantity

- 200 per package

Variations

Model No.	Description
FS083	Without Absorbing Sheet
FS084	With Absorbing Sheet

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI