

# Casting Material Mixing Bags

Durable mixing bags designed for preparing and handling casting materials in the field or lab. Made from tough, flexible plastic, they allow for easy mixing, pouring, and storage while minimizing mess and material loss. Ideal for forensic casting applications such as footprint or tire impression evidence.



## SPECIFICATIONS

### Construction

- Resealable
- 4 mm thick

### Dimensions

- 13"x18"

### Performance

- Allows a clean and simple way to mix water with casting material

### Quantity

- 100 per package

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI