


# Multi-Barrier Pouch

The Multi-barrier Pouch is a tamper-evident, seven-layer laminated pouch designed to protect FTA® cards from exposure to gases and liquids. It includes an outer paper layer for labeling and comes in small and large sizes.



## SPECIFICATIONS

### Construction

- Tamper-evident seal

### Quantity

- 100 per package

### Performance

- Outer paper layer can be used for labeling or writing
- Maintains sample security
- Protect the card from exposure to any gaseous or liquid external element

### Materials

- (7) laminated layers

## Variations

Model No.	Dimensions
FS161	3.7" X 3"
FS162	4.4" X 6.5"

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI