

Plastic-Polystyrene Jars, Various Sizes

Plastic polystyrene jars are lightweight, durable containers designed for secure collection, storage, and transport of forensic science samples. Available in various sizes, they provide a reliable, contamination-resistant option for preserving evidence such as trace materials, biological samples, and powders.



SPECIFICATIONS

Performance

- Ideal for fragile evidence such as paint chips and glass shards

Materials

- Plastic
- Polystyrene

Quantity

- 15 per case

Variations

| Model No. | Capacity |
|-----------|----------|
| FS163 | 1 oz |
| FS164 | 2 oz |
| FS165 | 4 oz |
| FS166 | 8 oz |
| FS167 | 16 oz |
| FS168 | 24 oz |
| FS169 | 32 oz |

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI