


# Bio-Foam® Impression Foam

Bio-Foam® Impression Foam is a soft, high-resolution casting material used in forensic science to capture detailed three-dimensional impressions of footprints, footwear, and other surface evidence. Its uniform compressibility allows it to conform precisely to fine surface features, preserving ridge detail and tread patterns for analysis.



FS170



FS171

## SPECIFICATIONS

### Construction

- Easy close lid
- Built-in storage for impressions

### Dimensions

- 12" L X 6" W X 2.25" H

### Performance

- Easy, clean, and efficient way to obtain an impression of suspected footwear

### Materials

- Sturdy cardboard
- Foam

## Variations

Model No.	Description
FS170	Single (1 Shoe)
FS171	Double (2 Shoes)

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI