

Spray Bottle With Adjustable Nozzle

A spray bottle with an adjustable nozzle is a versatile tool used in forensic investigations to apply liquids with controlled precision. The nozzle settings allow investigators to switch between a fine mist and a targeted stream, making it suitable for evenly distributing reagents such as luminol or other detection solutions.



FS172

SPECIFICATIONS

Construction

- Pump
- Adjustable nozzle

Capacity

- 16 oz

Performance

- Mist questioned areas

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI