

Shoe Covers

Lightweight yet durable, these shoe covers slip on easily to help prevent contamination from trace particles. Made from a woven, semi-impervious material, they feature an elastic top for a secure fit and anti-skid soles for traction on wet surfaces. Disposable and convenient, they're ideal for maintaining clean environments.



SPECIFICATIONS

Construction

- Elastic top band
- One size fits all
- Easily slips on over shoes
- Standard or Anti-Skid options available

Materials

- Semi-impervious, disposable Tyvek

Quantity

- 10 pairs

Variations

Model No.	Size
FS202	Anti-Skid
FS203	Standard

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI