

TranZport Hood

The TranZport™ Hood is a breathable, bacteria-filtering hood that limits contaminant spread while maintaining visibility and airflow. It has an elastic neckband and lock tab to prevent removal during transport.



FS213

SPECIFICATIONS

Construction

- Secure lock tab
- Elastic neckband
- Mesh netting

Quantity

- 5 per package

Performance

- Safely transport unruly or ill persons
- Fabric filters bacteria

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI