

# Bluestar® Magnum Formula Bloodstain Reagent

Invisible bloodstains produce an intense blue chemiluminescence (420 to 440 nanometer range) immediately after being sprayed with this product. Compared to traditional luminol reagent, this is three times brighter, is luminescent longer and is visible to the naked eye in semidarkness.



FS226

## SPECIFICATIONS

### Performance

- Ideal for miscellaneous items, surface, or clothes that have been washed
- Mix in a spray bottle containing 125 ml of distilled water
- Does not adversely affect DNA
- Odorless
- Fully soluble
- Makes invisible bloodstains produce an intense blue chemiluminescence

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI