

Bluestar Training Formula Bloodstain Reagent

Bluestar® Forensic Training is a practice version of the highly sensitive bloodstain detection reagent used in crime scene investigation. It produces the same visible blue chemiluminescent reaction for learning purposes, but it destroys DNA and is not suitable for forensic analysis.



SPECIFICATIONS

Performance

- For training purposes only
- Mix in a spray bottle containing 125 ml of distilled water
- Will destroy DNA and prevent DNA typing of the bloodstain
- Odorless
- Fully soluble
- Makes invisible bloodstains produce an intense blue chemiluminescence

Variations

Model No.	Quantity
FS227	4 per package
FS228	8 per package

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI