


EasiCollect Devices

These EasiCollect Devices are designed for consistent collection and transfer of cells onto included FTA® Cards using a foam collection surface. They securely hold the card in place with pressure clips to ensure uniform sample transfer and reliable results.



FS239

SPECIFICATIONS

Construction

- Collection surface
- Pressure clips

Materials

- Foam

Performance

- Compatible with manual or automated punching systems
- Captures a large number of cells for analysis

Quantity

- 50 per package

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI