

Reversible Backing Cards

These 3" x 5" Backing Cards are used to support lifted prints for documentation and analysis. Made from durable cast-coated cardstock with glossy and matte reversible sides, they allow both display and written labeling.



SPECIFICATIONS

Construction

- High gloss finish with matte reverse size

Dimensions

- 3" X 5"
- 0.012" thick

Materials

- Cast coated

Quantity

- 100 per pack

Variations

Model No.	Color
FS243	Black
FS244	White

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI