

# FTA Card

Designed for the collection, transportation, storage, and purification of nucleic acids from whole blood and buccal cells/saliva, each FTA® Card has an integral cover sheet which protects the sample and permits the card to stand on its own while drying. They come in multiple sizes with pre-defined sample areas.



## SPECIFICATIONS

### Performance

- Each circle holds 100–125  $\mu$ L
- For collection, transportation, storage, and purification of nucleic acids from blood and buccal cells/saliva

### Construction

- Integral cover sheet

### Quantity

- 100 per pack

## Variations

Model No.	Dimensions	Matrix	Circles
FS251	1.5" X 3"	1.5" X 1.5"	1
FS252	1.88" X 2.88"	1.25" X 2.88"	2
FS253	3" X 5"	3" X 3"	4

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI