

# FTA Elute Micro or Indicating Card

The FTA® Elute Card provides stable, long-term room temperature DNA storage with easy release using a simple elution method. It eliminates lengthy and multiple step isolation procedures, and the indicator card changes from purple to white when the card comes in contact with the sample.



## SPECIFICATIONS

### Construction

- Easy release using a simple elution method
- Integrated cover sheet

### Dimensions

- 1.5" X 3" with 1.5" X 1.5" matrix
- FS257: (4) 11 mm diameter circles
- FS258: (1) 1" diameter circle

### Quantity

- 100 per pack

### Performance

- Each circle holds 40 µL
- Indicator card changes from purple to white when it comes in contact with sample
- Eliminates lengthy and multiple step isolation procedures
- Provides stable, long-term room temperature DNA storage

## Variations

Model No.	Card Type	Circles
FS257	Elute	(4) 11 mm circles
FS258	Indicating Elute	(1) 1" circle

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI