

# Protein Saver 903 Card

The Protein Saver 903 Cards are a reliable sample collection and storage solution designed to preserve protein integrity during transport and long-term storage. Made with high-quality absorbent paper, it enables easy application and rapid drying of biological samples while maintaining stability at ambient temperatures.



FS259

## SPECIFICATIONS

### Construction

- (5) 0.5" circles

### Dimensions

- 3.5" X 2.67"

### Quantity

- 100 per package

### Performance

- Each circle holds 75 to 80  $\mu$ l
- Printed front and back
- Can be used with any type of fingerprinting ink

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI