

# Record Backing Card Pad

These Record Latent Print Backing Card Pads are used to mount and document single latent prints. Printed on glossy cardstock with designated print areas, they provide clear contrast for different powders.



## SPECIFICATIONS

### Construction

- High gloss print

### Materials

- Cardstock

### Dimensions

- Record card latent print space:  
(Smaller Card) 1.375" X 1.75"  
(Larger Card) 2.75" X 2.25"
- 0.012" thick

### Quantity

- 50/pad

## Variations

Model No.	Dimensions	Color
FS280	3" X 5"	Black
FS281	5.5" X 8.5"	Black
FS282	3" X 5"	White
FS283	5.5" X 8.5"	White

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI