

Epoxy Sample Containers

Epoxy Sample Containers are designed to preserve solid evidence with an airtight seal that prevents evaporation of flammables and hydrocarbons. The epoxy-coated interior resists rust, making them suitable for wet or damp items.



SPECIFICATIONS

Construction

- Prevents flammables and hydrocarbons from evaporating
- Airtight seal

Materials

- Coated with epoxy paint

Quantity

- 10 per case

Variations

Model No.	Capacity
FS321	0.5 Pint
FS322	1 Pint
FS323	1 Quart
FS324	1 Gallon

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI