

RIDGE Fiberglass Latent Print Powder Brush

The RIDGE™ Fiberglass Latent Print Powder Brush is used to apply fingerprint powder for developing latent prints. It has soft fiberglass bristles, a balanced plastic handle, and comes stored in an anti-static tube.



FS348

SPECIFICATIONS

Construction

- Print powder brush
- Anti-static tube

Dimensions

- Handle: 4.13" L
- Bristles: 2.25" L
- Tube: 7.5" L

Materials

- Fiberglass
- Plastic

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI