

Latent Print Powder Feather Duster

A high-quality marabou feather duster designed for the gentle and efficient application of latent fingerprint powders over large surfaces. Its soft, lightweight feathers allow controlled, even powder distribution on areas such as walls, doors, and countertops, helping preserve fine ridge detail while minimizing excess buildup.



SPECIFICATIONS

Dimensions

- Feathers: 5" L
- Overall Length: 9" L

Materials

- Marabou feathers

Quantity

- (1) Latent Print Powder Duster included

Performance

- For counters, walls, and doors

Variations

Model No.	Description
FS352	Black Duster
FS353	Red Duster
FS354	White Duster

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI