

Covered Dissection Table with Lid Assist, Elevating, Optional Ventilation

The Mopec Covered Dissection Tables combine safety, functionality, and discretion to meet the needs of anatomy labs and teaching facilities. Featuring the exclusive and secure Mopec lid assist, the hinged covers lock neatly under the table during use, minimizing disruptions and enhancing safety. Elevation ensures ergonomic comfort, making this table an ideal choice for fostering a safe and efficient dissection experience.



HB1500



HB1510



HB1550

SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- Built-in tapered slope
- Includes drain hole with PVC drain valve
- HB1500 and HB1510: (4) heavy duty 5" casters
- HB1550: (4) heavy duty 8" casters
- Casters consist of (3) total lock and (1) directional lock
- (4) Post sockets included

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 92" L x 30" W
- HB1500 and HB1550: 32"-39.9" H (hand crank elevation)
- HB1510: 35"-47" H
- HB1510: 6" Duct

Performance

- Mopec exclusive lid assist: hinged cover with swing latches for easy open and close
- Exhaust requirements: 220-340 CFM (HB1510)
- Static pressure: 0.34-0.67" WG (HB1510)
- ADA wheelchair accessible (HB1550)

Capacity

- 350 lbs

Options & Accessories

- Ventilation: Use HB1510
- Non-Elevation: See HB1400 series
- AX107: Bucket holder (Not compatible with ADA option)
- HO009: ADA Bucket Holder
- HO001: Bridge-style leg holder (must remove grid plate to work)
- HO005: Post leg-holders
- HO006: Post book holder
- HO024: Full-length shelf
- HO030: Instrument tray holder
- HO032: Short shelf (use with AX107)
- HO1110: Misting system

Related/Alternate/Obsolete Product Model Numbers

| Model No. | Description | Status | Notes |
|--------------|--|-----------|------------|
| 600021-HC-LA | Dissection Cart, Elevated, Lid Assist | Replaced | Use HB1500 |
| 600024-HC-LA | Dissection Cart, Ventilated, Elevated, Lid Assist | Replaced | Use HB1510 |
| HB1520 | Dissection Cart, Dual Ventilated, Elevated, Lid Assist | Obsoleted | Use HB1510 |

Variations

| Model No. | Description | Ventilation | Elevation |
|-----------|---------------------------------|-------------------------|----------------|
| HB1500 | Dissection Cart with Lid Assist | Non-Ventilated | Elevating |
| HB1510 | Dissection Cart with Lid Assist | Single Duct Ventilation | Elevating |
| HB1550 | Dissection Cart with Lid Assist | Non-Ventilated | Elevating, ADA |

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI