

Hydraulic Downdraft Dissection Table

Mopec's Hydraulic Downdraft Dissection Table provides a sturdy, easy-to-use surface for dissection work. It features downdraft ventilation to help control odors and airborne contaminants, along with hydraulic height adjustment for user comfort. Built with durable materials, it's well suited for labs, schools, and forensic settings.



HB300

SPECIFICATIONS

Construction

- Exhaust requirements: 350–500 CFM
- (4) heavy duty 8" casters
- Raise or lower the cart using the dual sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- The vented stainless-steel table allows fumes to be pulled down and away from the worksurface

Materials

- Vented 304 stainless steel table with No.4 polish and drainage outlet
- Hydraulic Stryker base

Dimensions

- 86" L (89" including duct) X 30" W X 30.25" to 43.75" H

Capacity

- 500 lbs

Utilities / Power

- Recommended airflow, adjustable between 350 and 500 CFM.
- Duct stub connection: 6" round
- Ball valve with 1" flexible drain

Options and Accessories

- HO001: Leg Rest Bridge

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI