

# Full Roller Pallet Assembly, Various Sizes

Mopec's Full Roller Pallet Assembly is engineered to provide smooth, controlled lifting and effortless transfer of heavy loads. For hydraulic forklifts, its full-roller surface ensures safe, efficient movement of trays, racks, or remains. Built with durable materials and designed for ergonomic handling.



M601H-23



M601H-31



M601H-27

## SPECIFICATIONS

### Construction

- M601H-23 and M601H-27: For use with corresponding JC Series trays (sold separately)
- M601H-31: For use with GA031 trays (sold separately)
- Galvanized steel rollers to allow for smooth movement of trays on the pallet
- Safety Latches with Bumpers on either side of the assembly to secure the tray in place during moving

### Materials

- 11ga 1.5" SQ 204 stainless steel tubing construction with a 180-grit finish

### Note:

- **Trays sold separately.**

### Variations

Model No.	Description	Dimensions
M601H-23	Full Roller Pallet Assembly for 23" Trays	81.8" L X 26.5" W X 6" H
M601H-27	Full Roller Pallet Assembly for 27" Trays	81.3" L X 31" W X 6" H
M601H-31	Full Roller Pallet Assembly for 31" Trays	81.3" L X 34" W X 6" H

### Dimensions

- Fork Pockets: 2.0" T X 4.0" W X 0.125" thick
- M601H-23: Assembly: 81.8" L X 26.5" W X 6.0" H  
Rollers: 1.4" diameter X 23.5" length
- M601H-27: Assembly: 81.3" L X 31" W X 6.0" H  
Rollers: 1.9" diameter X 28" length
- M601H-31: Assembly: 81.3" L X 34" W X 1.5" H  
Rollers: 1.9" diameter X 31" length

### Capacity

- 500 lbs

### Compatibility

- M601H-23: JC023 and JC023D trays
- M601H-27: JC027 and JC027D trays
- M601H-31: GA031 trays
- All compatible with M678 Hydraulic Lift JD405, JD450, and JD451 Forklifts with JD460 Forks
- Note: Not for use with hand crank lifts, use M600C Series instead.

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI