


Spray Hose Assembly with Backflow Protection

This spray hose assembly with backflow protection offers a convenient and reliable solution for spray cleaning in laboratory and clinical settings. Featuring an easy-grip handle with a flexible hose, it allows precise, comfortable operation. Designed to operate safely, it ensures effective cleaning and protects the water supply from contamination, making it ideal for maintaining a safe and efficient workspace.



OA-150

SPECIFICATIONS

Construction

- Easy grip assembly with flexible hose, conveniently placed for easy spray cleaning.
- Push button spray valve has 2.2 GPM aerator and 3/4-14 male thread.
- Maximum Pressure: 125 psi / Minimum Pressure: 20 psi

Material

- Flexible stainless steel hose assembly, 72" long.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI