

# Self-Propelled Covered Cadaver Carrier

Mopec's Self-Propelled Covered Cadaver Carrier is engineered to provide stable and dependable transport for heavy-duty mortuary applications. The hydraulic system allows smooth, controlled height adjustment. Designed for reliability and ease of use, it supports consistent performance in mortuary and pathology settings.



600039-SP

## SPECIFICATIONS

### Construction

- Includes non-removable tray top welded to cart
- Raise or lower the cart using the dual sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- Self-propelled drive system
- Removable and adjustable canopy assembly

### Materials

- 304 Stainless steel carrier top with a No. 4 polish finish
- Hydraulic Stryker base
- Canopy frame is made of rounded, lightweight aluminum
- Vinyl cover

### Capacity

- 1,000 lbs

### Dimensions

- 86" L X 38" W X 23.3" to 34.5" H (11" of elevation)
- Canopy Frame: 80" L X 28" W X 14.5" H
- Cover: 87" L X 41.5" W X 21" D

### Options & Accessories

- JZ030: Black vinyl cover (included and available individually)

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI