

Collection and Transport Bottles

These 100 ml Urine Specimen Collection Bottles feature a patented design that helps prevent tampering by deterring the use of substituted or externally heated specimens. Each bottle includes a concealed temperature monitoring thermometer hidden beneath a tear-strip label, allowing collectors to verify specimen temperature only after collection.



FS005

SPECIFICATIONS

Construction

- Concealed temperature monitoring thermometer
- Heat-sealed integrity pouch
- External label w/ tear-strip

Performance

- Transport and collect urine samples
- Prevents subjects from using someone else's urine

Capacity

- 100 ml

Quantity

- 100 per case

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI