

Heavy-Duty Body Bags

These heavy-duty body bags are made from durable, chlorine-free polyethylene that remains flexible in extreme temperatures and supports substantial weight. Designed with multiple reinforced carry straps and a secure inner liner, they provide reliable handling and effective containment while meeting OSHA standards for safety and protection.



FS074

SPECIFICATIONS

Construction

- Equipped with 10 carry straps
- Heat-sealed inner bag liner

Weight Capacity

- 450 lbs

Performance

- Complies with OSHA regulation 3130
- Long shelf life
- Stays flexible up to - 70° F

Quantity

- FS074: 5 per case
- FS075: 3 per case

Materials

- 12 mm chlorine-free polyethylene

Variations

Model No.	Dimensions	Quantity
FS074	36" X 96"	5 per case
FS075	48" X 100"	3 per case

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI