


Impression Foam

Impression Foam is a soft, high-resolution casting material used in forensic science to capture detailed three-dimensional impressions of footprints, footwear, and other surface evidence. Its uniform compressibility allows it to conform precisely to fine surface features, preserving ridge detail and tread patterns for analysis.



FS170



FS171

SPECIFICATIONS

Construction

- Easy close lid
- Built-in storage for impressions

Dimensions

- 12" L X 6" W X 2.25" H

Performance

- Easy, clean, and efficient way to obtain an impression of suspected footwear

Materials

- Sturdy cardboard
- Foam

Variations

Model No.	Description
FS170	Single (1 Shoe)
FS171	Double (2 Shoes)

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI