

Crime Scene Eye Forensic Glasses

Lightweight anti-fog polycarbonate goggles provide 100% UV protection and reduce eye fatigue and blue haze. They meet ANSI Z87.1 and OSHA 1910.133 standards and fit over prescription glasses. Available in clear, orange, yellow, and red lenses.



SPECIFICATIONS

Construction

- Interchangeable lenses
- Non-slip nose grip

Quantity

- 2 per pack

Materials

- Polycarbonate

Performance

- Allows fluorescence to be seen

Variations

Model No.	Description	Block Wavelengths
FS179	Orange Eye Glasses	Below 570nm
FS180	Red Eye Glasses	Below 610nm
FS181	Yellow Eye Glasses	Below 500nm
FS182	Yellow/Orange Eye Glasses	Below 500nm/570nm

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI