

# Mini Forensic Kit

The Mini Forensic Kit includes all the essentials needed to detect bloodstains at a crime scene. The tablets dissolve in the pre-filled spray bottle within 2–3 minutes, making the reagent quickly ready for use. It offers a simple, compact solution for reliable field application.



FS225

## SPECIFICATIONS

### Construction

- Kit consists of:
  - (2) Pre-filled spray bottles
  - (2) Foil pouches
- Each foil pouch has (2) tablets

### Performance

- Everything needed for locating bloodstains at a crime scene
- Drop the tablets into one of the pre-filled spray bottles
- Tablets will dissolve within 2–3 minutes and will be ready to use

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI