

Bluestar Mini Kit

The Bluestar® Mini Kit includes all the essentials needed to detect bloodstains at a crime scene. The tablets dissolve in the pre-filled spray bottle within 2–3 minutes, making the reagent quickly ready for use. It offers a simple, compact solution for reliable field application.



FS225

SPECIFICATIONS

Construction

- Kit consists of:
 - (2) Pre-filled spray bottles
 - (2) Foil pouches
- Each foil pouch has 2 Bluestar tablets

Performance

- Everything needed for locating bloodstains at a crime scene
- Drop the Bluestar tablets into one of the pre-filled spray bottles
- Bluestar tablets will dissolve within 2–3 minutes and will be read to use

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI