

FTA Indicating Card

FTA® Indicating Cards are used to collect and store buccal cell and saliva samples for DNA analysis, changing from pink to white when sample is applied. They include an integral cover sheet for protection and stability during drying, with multiple size options for different sample capacities.



SPECIFICATIONS

Performance

- Color changes from pink to white when card comes into contact with samples
- For collection and storage of buccal cell and saliva samples for DNA analysis
- Each circle holds 100–125 µL

Construction

- Integral cover sheet

Quantity

- 100 per pack

Variations

Model No.	Dimensions	Matrix	Circles
FS254	1.5" X 3"	1.5" X 1.5"	1
FS255	1.88" X 2.88"	1.25" X 2.88"	2
FS256	3" X 5"	3" X 3"	4

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI