

# Solid Sample Metal Evidence Containers, Various Sizes

These Solid Sample Metal Evidence Containers are used for preserving solid forms of physical evidence. Constructed as metal containers with secure lids, they form an airtight seal once closed, preventing flammables and hydrocarbons from evaporating into the atmosphere. Designed as a standard method for evidence storage, they help maintain the integrity of collected samples.



## SPECIFICATIONS

### Construction

- Metal containers and lids of various sizes to hold evidence from crime scenes
- (3) size options available

### Quantity

- 10 per case

### Materials

- Metal

### Performance

- Collects and preserves “solid” forms of physical evidence
- Air tight seal is formed once the lid is on which prevents flammables and hydrocarbons from evaporating into the air

## Variations

Model No.	Description	Size
FS325	Solid Sample Metal Evidence Containers	1 Pint
FS326	Solid Sample Metal Evidence Containers	1 Quart
FS327	Solid Sample Metal Evidence Containers	1 Gallon

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI