

Sterile Fiberglass Latent Print Powder Brush

The Sterile Fiberglass Latent Print Powder Brush is used to apply fingerprint powder while maintaining a sterile, contamination-free condition. It features a plastic handle and fiberglass bristles, and is individually packaged for maximum safety and compliance with medical packaging standards.



FS349

SPECIFICATIONS

Construction

- Transparent packaging
- Print powder brush
- Conforms with EN-868-5
- Sterilization Technique: Ethylene Oxide 15 gr EtO ampoule is used for per 35 dm³ container

Dimensions

- Handle: 4.13" L
- Brush: 1.77" L

Materials

- Co-polymer film
- Craft paper
- Plastic

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI