

# Blood Tube Specimen Security Seal

These Blood Tube Specimen Security Seal are designed to safeguard sample integrity from collection to analysis. This tamper-evident seal securely wraps around blood tubes, providing a clear visual indication if a specimen has been opened or compromised.



FS359

## SPECIFICATIONS

### Performance

- Tamper evident, leaves a void print if there's an attempt to remove
- Can be easily written on

### Dimensions

- 7" L X 1" W

### Quantity

- 100 per pack

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI