

Slab and Roller Cleaner

Slab and Roller Cleaner is a powerful, fast-acting solution designed to remove inks, grease, and buildup from printing slabs and rollers. Its effective formula cuts through tough residues while helping maintain the performance and longevity of equipment.



SPECIFICATIONS

Construction

- Aerosol spray can
- Fast acting ink solvent

Performance

- Loosens and removes fingerprint ink from rubber rollers, glass or stainless steel slabs

Variations

Model No.	Description
FS367	Cleaner, 6 oz
FS368	Cleaner, 10 oz

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI