

Casting Frames

Casting Frames are durable tools designed to securely hold impression materials in place during the casting process. They help maintain proper shape and stability while materials set, ensuring accurate and consistent results. Tire print or Footprint options available.



SPECIFICATIONS

Materials

- 0.04" thick aluminum

Dimensions

- FS388: 12-20" L X 8.5" X 1.75" H
- FS389: 20.25-24" L X 12" W X 1.75" H

Performance

- Adjustable
- Restricts the flow and controls the depth of casting material

Variations

Model No.	Description	Dimensions
FS388	Footprint Casting Frame	12-20" L X 8.5" X 1.75" H
FS389	Tire Print Casting Frame	20.25-24" L X 12" W X 1.75" H

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI